

TRUE METRIX® PRO

The McKesson TRUE METRIX® PRO Professional Blood Glucose Monitoring System is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertip or forearm and venous whole blood. It is intended and approved for multiple-patient use in professional healthcare settings such as long term care facilities and physician offices. The meter features Triple Sense Technology™ and offers enhanced lifestyle management features, such as event tags and testing reminders.

Features and Benefits

- No coding required
- Results as fast as 4 seconds
- Tiny 0.5 microliter sample size
- Alternate site testing (fingertip or forearm)
- Stores up to 500 results with time/date
- Event tagging; analyzes environmental and physiological factors to better interpret tests
- 4 testing reminder alarms
- Audible fill detection for blood or control sample
- Ketone test reminder
- Easy-to-use "hands free" strip release button



ORDERING INFORMATION

MFR #	DESCRIPTION	QTY
06-RE4051P-00	TRUE METRIX® PRO Monitoring System Includes: Meter, Installed 3-Volt Battery, Owner's Manual, Self-Test Log Book, Carry Case	1 ea, 6 ea/cs
06-R3051P-01	TRUE METRIX® PRO Strips, 100 count	50 ea/vial, 2 vials/bx, 12 bx/cs
06-R3051P-05	TRUE METRIX® PRO Strips, 50 count	50 ea/vial, 1 vial/bx, 24 bx/cs
06-R5051-1	TRUE METRIX® Control Level 1, 70-90 mg/dL	1 ea/bx, 24 bx/cs
06-R5051-2	TRUE METRIX® Control Level 2, 150-200 mg/dL	1 ea/bx, 24 bx/cs
06-R5051-3	TRUE METRIX® Control Level 3, 301-398 mg/dL	1 ea/bx, 24 bx/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.